





RDX Antibody

Product Code	CSB-PA120056
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35241
Immunogen	Fusion protein of Human RDX
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	radixin
Immunogen Species	Homo sapiens (Human)
Target Names	RDX


Image



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA120056(RDX Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA120056(RDX Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 6%SDS-PAGE, Lysate: 40 μ g, Lane 1-5:
Hela cells, human brain tissue, 293T cells,
human hepatocellular carcinoma tissue, NIH/3T3
cells, Primary antibody: CSB-PA120056(RDX
Antibody) at dilution 1/600, Secondary antibody:
Goat anti rabbit IgG at 1/8000 dilution, Exposure
time: 5 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.