





FSHR Antibody

Product Code	CSB-PA288553
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23945
Immunogen	Fusion protein of Human FSHR
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:2000-1:5000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	follicle stimulating hormone receptor
Immunogen Species	Homo sapiens (Human)
Target Names	FSHR

Image



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA288553(FSHR Antibody) at dilution 1/70, on the right is treated with fusion protein. (Original magnification: ×200)



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA288553(FSHR Antibody) at dilution 1/70, on the right is treated with fusion protein. (Original magnification: ×200)