





GNA11 Antibody

Product Code	CSB-PA201320
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29992
Immunogen	Fusion protein of Human GNA11
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	guanine nucleotide binding protein (G protein), alpha 11 (Gq class)
Immunogen Species	Homo sapiens (Human)
Target Names	GNA11

Image



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA201320(GNA11 Antibody) at dilution 1/40, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane: Mouse eyes tissue, Primary antibody: CSB-PA201320(GNA11 Antibody) at dilution 1/550, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 10 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.