





GNAS Antibody

Product Code	CSB-PA247539
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63092
Immunogen	Fusion protein of Human GNAS
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:5000,WB:1:200-1:1000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	GNAS complex locus
Immunogen Species	Homo sapiens (Human)
Target Names	GNAS

Image



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA247539(GNAS Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)



Gel: 10%SDS-PAGE, Lysate: 30 μg, Lane: Mouse brain tissue, Primary antibody: CSB-PA247539(GNAS Antibody) at dilution 1/500, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 5 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.