




GADD45G Antibody

Product Code	CSB-PA620207
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95257
Immunogen	Fusion protein of Human GADD45G
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:2000-1:10000,IHC:1:100-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	growth arrest and DNA-damage-inducible, gamma
Immunogen Species	Homo sapiens (Human)
Target Names	GADD45G

Image



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA620207(GADD45G Antibody) at dilution 1/120, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.