




GDF15 Antibody

Product Code	CSB-PA255283
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99988
Immunogen	Fusion protein of Human GDF15
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:1000-1:2000,IHC:1:10-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	growth differentiation factor 15
Immunogen Species	Homo sapiens (Human)
Target Names	GDF15

Image



The image on the left is immunohistochemistry of paraffin-embedded Human renal cancer tissue using CSB-PA255283(GDF15 Antibody) at dilution 1/7, on the right is treated with fusion protein. (Original magnification: x200)