





MGEA5 Antibody

Product Code	CSB-PA910720
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60502
Immunogen	Fusion protein of Human MGEA5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:2000,WB:1:200-1:1000,IHC:1:15-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	Meningioma expressed antigen 5 (hyaluronidase)
Immunogen Species	Homo sapiens (Human)
Target Names	MGEA5

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA910720(MGEA5 Antibody) at dilution 1/15, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane: Mouse pancreas tissue, Primary antibody: CSB-PA910720(MGEA5 Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.