




HIST1H1C Antibody

Product Code	CSB-PA296044
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16403
Immunogen	Fusion protein of Human HIST1H1C
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:2000-1:10000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	histone cluster 1, H1c
Immunogen Species	Homo sapiens (Human)
Target Names	HIST1H1C

Image



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA296044(HIST1H1C Antibody) at dilution 1/70, on the right is treated with fusion protein. (Original magnification: x200)