




HSPA1A Antibody

Product Code	CSB-PA137102
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DMV8
Immunogen	Fusion protein of Human HSPA1A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:2000-1:5000,IHC:1:50-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	heat shock 70kDa protein 1A
Immunogen Species	Homo sapiens (Human)
Target Names	HSPA1A

Image



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA137102(HSPA1A Antibody) at dilution 1/40, on the right is treated with fusion protein. (Original magnification: x200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.