




# IGF2BP3 Antibody

<b>Product Code</b>	CSB-PA745585
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O00425
<b>Immunogen</b>	Fusion protein of Human IGF2BP3
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse
<b>Tested Applications</b>	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Alias</b>	insulin-like growth factor 2 mRNA binding protein 3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	IGF2BP3

## Image



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane 1-3: A172 cells, hela cells, hepg2 cells, Primary antibody: CSB-PA745585(IGF2BP3 Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 20 seconds

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.