





INHBA Antibody

Product Code	CSB-PA903855
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08476
Immunogen	Fusion protein of Human INHBA
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	inhibin, beta A
Immunogen Species	Homo sapiens (Human)
Target Names	INHBA


Image



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA903855(INHBA Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA903855(INHBA Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 8%SDS-PAGE, Lysate: 40 μ g, Lane: Mouse fat tissue, Primary antibody: CSB-PA903855(INHBA Antibody) at dilution 1/701, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute