





ITGB3BP Antibody

Product Code	CSB-PA980626
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13352
Immunogen	Fusion protein of Human ITGB3BP
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC;ELISA:1:1000-1:5000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	integrin beta 3 binding protein (beta3-endonexin)
Immunogen Species	Homo sapiens (Human)
Target Names	ITGB3BP

Image



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using CSB-PA980626(ITGB3BP Antibody) at dilution 1/50. (Original magnification: ×200)



The image is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA980626(ITGB3BP Antibody) at dilution 1/50. (Original magnification: ×200)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.