





IRF9 Antibody

Product Code	CSB-PA570795
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00978
Immunogen	Fusion protein of Human IRF9
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; ELISA: 1:1000-1:2000, WB: 1:200-1:1000, IHC: 1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	interferon regulatory factor 9
Immunogen Species	Homo sapiens (Human)
Target Names	IRF9


Image



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA570795 (IRF9 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using CSB-PA570795 (IRF9 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 8%SDS-PAGE, Lysate: 40 μ g, Lane 1-2:
Mouse panceas tissue, NIH/3T3 cells, Primary
antibody: CSB-PA570795(IRF9 Antibody) at
dilution 1/380, Secondary antibody: Goat anti
rabbit IgG at 1/8000 dilution, Exposure time: 1.5
minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.