




MAPK8 Antibody

Product Code	CSB-PA236755
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45983
Immunogen	Fusion protein of Human MAPK8
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	mitogen-activated protein kinase 8
Immunogen Species	Homo sapiens (Human)
Target Names	MAPK8

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane: A172 cells, Primary antibody: CSB-PA236755(MAPK8 Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 3 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.