




KCNA5 Antibody

Product Code	CSB-PA229768
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22460
Immunogen	Fusion protein of Human KCNA5
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	potassium voltage-gated channel, shaker-related subfamily, member 5
Immunogen Species	Homo sapiens (Human)
Target Names	KCNA5

Image



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane: Mouse spleen tissue, Primary antibody: CSB-PA229768(KCNA5 Antibody) at dilution 1/700, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 60 seconds