




BIRC7 Antibody

Product Code	CSB-PA829257
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CA5
Immunogen	Full length fusion protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	baculoviral IAP repeat containing 7
Immunogen Species	Homo sapiens (Human)
Target Names	BIRC7

Image



Gel: 8%SDS-PAGE,Lysate: 40 µg,Lane 1-4: A549 cells, Jurkat cells, Raji cells, HeLa cells,Primary antibody: CSB-PA829257(BIRC7 Antibody) at dilution 1/600 dilution,Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,Exposure time: 25 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.