





MAGED1 Antibody

Product Code	CSB-PA695195
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5V3
Immunogen	Fusion protein of Human MAGED1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:5000,WB:1:1000-1:5000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	melanoma antigen family D, 1
Immunogen Species	Homo sapiens (Human)
Target Names	MAGED1


Image



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA695195(MAGED1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA695195(MAGED1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)



Gel: 10%SDS-PAGE, Lysate: 40 μ g, Lane:
Human fetal brain tissue, Primary antibody:
CSB-PA695195(MAGED1 Antibody) at dilution
1/1200, Secondary antibody: Goat anti rabbit IgG
at 1/8000 dilution, Exposure time: 90 seconds