





# MAGED1 Antibody

<b>Product Code</b>	CSB-PA695195
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9Y5V3
<b>Immunogen</b>	Fusion protein of Human MAGED1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB,IHC;ELISA:1:1000-1:5000,WB:1:1000-1:5000,IHC:1:25-1:100
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Alias</b>	melanoma antigen family D, 1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MAGED1


## Image



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA695195(MAGED1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA695195(MAGED1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 40  $\mu$ g, Lane:  
Human fetal brain tissue, Primary antibody:  
CSB-PA695195(MAGED1 Antibody) at dilution  
1/1200, Secondary antibody: Goat anti rabbit IgG  
at 1/8000 dilution, Exposure time: 90 seconds

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.