




TMPRSS7 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA057341 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q7RTY8 |
| Immunogen | Fusion protein of Human TMPRSS7 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | -20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |
| Purification Method | Antigen affinity purification |
| Isotype | IgG |
| Alias | transmembrane protease, serine 7 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | TMPRSS7 |

Image



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane 1-4: A549 cells, hepg2 cells, lovo cells, A431 cells, Primary antibody: CSB-PA057341(TMPrSS7 Antibody) at dilution 1/250, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 5 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.