




MFN1 Antibody

Product Code	CSB-PA929799
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWA4
Immunogen	Fusion protein of Human MFN1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:2000-1:5000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	mitofusin 1
Immunogen Species	Homo sapiens (Human)
Target Names	MFN1

Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA929799(MFN1 Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: x200)