





RPSA Antibody

Product Code	CSB-PA937841
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08865
Immunogen	Fusion protein of Human RPSA
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	ribosomal protein SA
Immunogen Species	Homo sapiens (Human)
Target Names	RPSA


Image



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA937841(RPSA Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA937841(RPSA Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 40 μ g, Lane 1-5:
A549 cells, NIH/3T3 cells, 293T cells, human
hepatocellular carcinoma tissue, hela cells,
Primary antibody: CSB-PA937841(RPSA
Antibody) at dilution 1/433.3, Secondary
antibody: Goat anti rabbit IgG at 1/8000 dilution,
Exposure time: 15 seconds