





MSH6 Antibody

Product Code	CSB-PA051370
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52701
Immunogen	Fusion protein of Human MSH6
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:2000,WB:1:200-1:1000,IHC:1:10-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	mutS homolog 6 (E. coli)
Immunogen Species	Homo sapiens (Human)
Target Names	MSH6

Image



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA051370(MSH6 Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-2: HeLa cells, 293T cells, Primary antibody: CSB-PA051370(MSH6 Antibody) at dilution 1/300, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 45 seconds