





MYH7 Antibody

Product Code	CSB-PA665912
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12883
Immunogen	Fusion protein of Human MYH7
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,IHC;ELISA:1:1000-1:5000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	myosin, heavy chain 7, cardiac muscle, beta
Immunogen Species	Homo sapiens (Human)
Target Names	MYH7

Image



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using CSB-PA665912(MYH7 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA665912(MYH7 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)