




NME3 Antibody

Product Code	CSB-PA058575
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13232
Immunogen	Fusion protein of Human NME3
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;ELISA:1:1000-1:5000,WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	NME/NM23 nucleoside diphosphate kinase 3
Immunogen Species	Homo sapiens (Human)
Target Names	NME3

Image



Gel: 12%SDS-PAGE, Lysate: 40 µg, Lane: Mouse kidney tissue, Primary antibody: CSB-PA058575(NME3 Antibody) at dilution 1/800, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.