





NOS2 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA545107 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P35228 |
| Immunogen | Fusion protein of Human NOS2 |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA,IHC;ELISA:1:2000-1:10000,IHC:1:100-1:300 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | -20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |
| Purification Method | Antigen affinity purification |
| Isotype | IgG |
| Alias | nitric oxide synthase 2, inducible |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | NOS2 |

Image



The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using CSB-PA545107(NOS2 Antibody) at dilution 1/40, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA545107(NOS2 Antibody) at dilution 1/40, on the right is treated with fusion protein. (Original magnification: x200)