




BIRC2 Antibody

Product Code	CSB-PA191509
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13490
Immunogen	Fusion protein of Human BIRC2
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; ELISA: 1:1000-1:2000, WB: 1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	BIRC2

Image



Gel: 8% SDS-PAGE, Lysate: 40 µg, Lane: Human fetal muscle tissue, Primary antibody: CSB-PA191509 (BIRC2 Antibody) at dilution 1/300, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 2 minutes

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.