




PGBD2 Antibody

Product Code	CSB-PA214066
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P3X8
Immunogen	Fusion protein of Human PGBD2
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	PGBD2

Image



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane 1-4: A375 cells, A172 cells, mouse pancreas tissue, human liver cancer tissue, Primary antibody: CSB-PA214066(PGBD2 Antibody) at dilution 1/550, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 20 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.