





AKAP7 Antibody

Product Code	CSB-PA097902
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43687
Immunogen	Fusion protein of Human AKAP7
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:2000,WB:1:200-1:1000,IHC:1:15-1:50
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	AKAP7


Image



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA097902(AKAP7 Antibody) at dilution 1/15, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA097902(AKAP7 Antibody) at dilution 1/15, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 12%SDS-PAGE, Lysate: 40 μ g, Lane:
 Human bladder carcinoma tissue, Primary
 antibody: CSB-PA097902(AKAP7 Antibody) at
 dilution 1/200 dilution, Secondary antibody:
 Goat anti rabbit IgG at 1/8000 dilution, Exposure
 time: 30 seconds