





AKAP8L Antibody

Product Code	CSB-PA547116
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULX6
Immunogen	Fusion protein of Human AKAP8L
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:200-1:1000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	AKAP8L


Image



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA547116(AKAP8L Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA547116(AKAP8L Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 8%SDS-PAGE, Lysate: 40 μ g, Lane 1-3:
Jurkat cells, 293T cells, K562 cells, Primary
antibody: CSB-PA547116(AKAP8L Antibody) at
dilution 1/117, Secondary antibody: Goat anti
rabbit IgG at 1/8000 dilution, Exposure time: 30
seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.