





ANKRD1 Antibody

Product Code	CSB-PA562564
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15327
Immunogen	Fusion protein of Human ANKRD1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	ANKRD1


Image



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA562564(ANKRD1 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using CSB-PA562564(ANKRD1 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 8%SDS-PAGE, Lysate: 60 μ g, Lane 1-3:
 Mouse lung tissue, NIH/3T3 cells, hepg2 cells,
 Primary antibody: CSB-PA562564(ANKRD1
 Antibody) at dilution 1/400, Secondary antibody:
 Goat anti rabbit IgG at 1/8000 dilution, Exposure
 time: 3 seconds