





IDH3B Antibody

Product Code	CSB-PA114986
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43837
Immunogen	Fusion protein of Human IDH3B
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:1000-1:2000,WB:1:200-1:1000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	IDH3B

Image



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA114986(IDH3B Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 40 µg, Lane 1-5: HepG2 cells, mouse kidney tissue, lovo cells, mouse eyes tissue, hela cells, Primary antibody: CSB-PA114986(IDH3B Antibody) at dilution 1/233.3, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 20 seconds