





SIGMAR1 Antibody

Product Code	CSB-PA233715
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99720
Immunogen	Fusion protein of Human SIGMAR1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	SIGMAR1


Image



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using CSB-PA233715(SIGMAR1 Antibody) at dilution 1/50. (Original magnification: x200)



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA233715(SIGMAR1 Antibody) at dilution 1/50. (Original magnification: x200)



Gel: 10%SDS-PAGE, Lysate: 40 μ g, Lane 1-3:
Mouse liver tissue, Mouse brain tissue, A549
cells, Primary antibody: CSB-
PA233715(SIGMAR1 Antibody) at dilution 1/310,
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution, Exposure time: 1 second

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.