





AMY2A Antibody

Product Code	CSB-PA201757
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04746
Immunogen	Fusion protein of Human AMY2A
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:10000,WB:1:1000-1:5000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	AMY2A


Image



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using CSB-PA201757(AMY2A Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using CSB-PA201757(AMY2A Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)



Gel: 15%SDS-PAGE, Lysate: 40 μ g, Lane:
Mouse pancreas tissue, Primary antibody: CSB-
PA201757(AMY2A Antibody) at dilution 1/1400,
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution, Exposure time: 1 second