




MPG Antibody

Product Code	CSB-PA832374
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29372
Immunogen	Fusion protein of Human MPG
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:2000-1:10000,WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	MPG

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-7: 231 cells, 293T cells, Lovo cells, hepG2 cells, A549 cells, PC3 cells, A172 cells, Primary antibody: CSB-PA832374(MPG Antibody) at dilution 1/1500, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.