




# CYP11A1 Antibody

<b>Product Code</b>	CSB-PA048582
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P05108
<b>Immunogen</b>	Fusion protein of Human CYP11A1
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; ELISA: 1:1000-1:2000, WB: 1:200-1:1000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
<b>Purification Method</b>	Antigen affinity purification
<b>Isotype</b>	IgG
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	CYP11A1

## Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane:  
Human normal kidney tissue, Primary antibody:  
CSB-PA048582(CYP11A1 Antibody) at dilution  
1/200, Secondary antibody: Goat anti rabbit IgG  
at 1/8000 dilution, Exposure time: 5 minutes