





FAM89B Antibody

Product Code	CSB-PA290518
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5H3
Immunogen	Fusion protein of Human FAM89B
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:200-1:1000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	FAM89B


Image



The image on the left is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using CSB-PA290518(FAM89B Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA290518(FAM89B Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)



Gel: 12%SDS-PAGE, Lysate: 40 μ g, Lane:
 Mouse liver tissue, Primary antibody: CSB-
 PA290518(FAM89B Antibody) at dilution 1/400,
 Secondary antibody: Goat anti rabbit IgG at
 1/8000 dilution, Exposure time: 10 seconds