





FRK Antibody

Product Code	CSB-PA570140
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42685
Immunogen	Fusion protein of Human FRK
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	FRK


Image



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA570140(FRK Antibody) at dilution 1/50. (Original magnification: ×200)



The image is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA570140(FRK Antibody) at dilution 1/50. (Original magnification: ×200)



Gel: 6%SDS-PAGE, Lysate: 40 μ g, Lane 1-2:
K562 cells, A431 cells, Primary antibody: CSB-
PA570140(FRK Antibody) at dilution 1/1200,
Secondary antibody: Goat anti rabbit IgG at
1/8000 dilution, Exposure time: 40 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.