




MDH1 Antibody

Product Code	CSB-PA449427
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40925
Immunogen	Fusion protein of Human MDH1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:2000-1:5000,WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	MDH1

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-3: Mouse heart tissue, Mouse liver tissue, Mouse stomach tissue, Primary antibody: CSB-PA449427(MDH1 Antibody) at dilution 1/500, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 3 seconds