




OVAL Antibody

Product Code	CSB-PA271972
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01012
Immunogen	OVAL protein
Raised In	Rabbit
Species Reactivity	Chicken
Tested Applications	ELISA, WB; ELISA: 1:2000-1:5000, WB: 1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Gallus gallus (Chicken)
Target Names	SERPINB14

Image



Gel: 10%SDS-PAGE, Lysate: 0.4 µg, 0.8 µg, Lane 1-2: 0.4 µg, 0.8 µg OVAL protein, Primary antibody: CSB-PA271972(OVAL Antibody) at dilution 1/500, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 3 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.