




MEP1A Antibody

Product Code	CSB-PA047131
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16819
Immunogen	Fusion protein of Human MEP1A
Raised In	Rabbit
Species Reactivity	Human,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	MEP1A

Image



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane: Human rectum carcinoma tissue, Primary antibody: CSB-PA047131(MEP1A Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 3 minutes

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.