





REG3G Antibody

Product Code	CSB-PA074697
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UW15
Immunogen	Fusion protein of Human REG3G
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:100-1:300
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	REG3G

Image



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using CSB-PA074697(REG3G Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: X200)



Gel: 10%SDS-PAGE, Lysate: 40μg, Lane: Mouse pancreas tissue, Primary antibody: CSB-PA074697(REG3G Antibody) at dilution 1/400, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 40 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.