




HMGN5 Antibody

Product Code	CSB-PA944763
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P82970
Immunogen	Fusion protein of Human HMGN5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	HMGN5

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-4:
 Raji cells, 293T cells, Hela cells, 231 cells,
 Primary antibody: CSB-PA944763(HMGN5
 Antibody) at dilution 1/200, Secondary antibody:
 Goat anti rabbit IgG at 1/8000 dilution, Exposure
 time: 5 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.