




NUMB Antibody

Product Code	CSB-PA909836
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49757
Immunogen	Fusion protein of Human NUMB
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	NUMB

Image



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane 1-6: K562 cells, hela cells, skov3 cells, NIH/3T3 cells, Raji cells, human testis tissue, Primary antibody: CSB-PA909836(NUMB Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 2 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.