




NR4A1 Antibody

Product Code	CSB-PA179677
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22736
Immunogen	Fusion protein of Human NR4A1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	nuclear receptor subfamily 4, group A, member 1
Immunogen Species	Homo sapiens (Human)
Target Names	NR4A1

Image



Gel: 6%SDS-PAGE, Lysate: 40 µg, Lane: PC3 cells, Primary antibody: CSB-PA179677(NR4A1 Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 5 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.