





DSTYK Antibody

Product Code	CSB-PA846981
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XUX3
Immunogen	Fusion protein of Human DSTYK
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; ELISA: 1:1000-1:2000, WB: 1:200-1:1000, IHC: 1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Alias	dual serine/threonine and tyrosine protein kinase
Immunogen Species	Homo sapiens (Human)
Target Names	DSTYK

Image



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA846981 (DSTYK Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 6% SDS-PAGE, Lysate: 40 µg, Lane 1-2: 823 cells, hepg2 cells, Primary antibody: CSB-PA846981 (DSTYK Antibody) at dilution 1/200, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 1 minute

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.