





RPS3 Antibody

Product Code	CSB-PA152933
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23396
Immunogen	Full length fusion protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	RPS3

Image



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA152933(RPS3 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)



Gel: 12%SDS-PAGE,Lysate: 40 µg,Lane 1-7: HepG2 cells, Jurkat cells, Hela cells, Raji cells, RAW264.7 cells, A549 cells, 293T cells,Primary antibody: CSB-PA152933(RPS3 Antibody) at dilution 1/350 dilution,Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,Exposure time: 10 seconds