




S100PBP Antibody

Product Code	CSB-PA125380
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BU1
Immunogen	Fusion protein of Human S100PBP
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB;ELISA:1:1000-1:2000,WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	S100PBP

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane:
Human fetal brain tissue, Primary antibody:
CSB-PA125380(S100PBP Antibody) at dilution
1/350, Secondary antibody: Goat anti rabbit IgG
at 1/8000 dilution, Exposure time: 1 minute