





STX10 Antibody

Product Code	CSB-PA225382
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60499
Immunogen	Fusion protein of Human STX10
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:10000,WB:1:1000-1:5000,IHC:1:30-1:150
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	STX10


Image



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using CSB-PA225382(STX10 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA225382(STX10 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 12%SDS-PAGE, Lysate: 40 μ g, Lane 1-4:
A375 cells, A431 cells, lovo cells, hela cells,
Primary antibody: CSB-PA225382(STX10
Antibody) at dilution 1/800, Secondary antibody:
Goat anti rabbit IgG at 1/8000 dilution, Exposure
time: 5 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.