





TGM5 Antibody

Product Code	CSB-PA965532
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43548
Immunogen	Fusion protein of Human TGM5
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:5000,WB:1:500-1:2000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	TGM5


Image



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using CSB-PA965532(TGM5 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using CSB-PA965532(TGM5 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 6%SDS-PAGE,Lysate: 40 μ g,Lane 1-2:
Mouse heart tissue, MCF7 cells,Primary
antibody: CSB-PA965532(TGM5 Antibody) at
dilution 1/400 dilution,Secondary antibody: Goat
anti rabbit IgG at 1/8000 dilution,Exposure time:
2 minutes

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.