




UQCRC2 Antibody

Product Code	CSB-PA569709
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22695
Immunogen	Fusion protein of Human UQCRC2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;ELISA:1:2000-1:10000,WB:1:1000-1:5000,IHC:1:25-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification Method	Antigen affinity purification
Isotype	IgG
Immunogen Species	Homo sapiens (Human)
Target Names	UQCRC2


Image



The image on the left is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using CSB-PA569709(UQCRC2 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using CSB-PA569709(UQCRC2 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x200)



Gel: 8%SDS-PAGE, Lysate: 40 μ g, Lane 1-5:
Jurkat, Raji, hepg2, A431 and HT-29 cell,
Primary antibody: CSB-PA569709(UQCRC2
Antibody) at dilution 1/400 dilution, Secondary
antibody: Goat anti rabbit IgG at 1/8000 dilution,
Exposure time: 10 seconds

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.